

## The Blocking Solution

(article no. 110)

## Casein based blocker for saturating free binding sites on plastic surfaces or other protein binding surfaces and for minimizing nonspecific binding to surfaces

	For research use only, not for diagnostic use
Expiry date when stored unopened:	please refer to the label on the bottle
Preservative:	contains < 0.0014 % [w/w] reaction mass of CMIT/MIT (3:1)
pH-value:	$7.2 \pm 0.2$
Storage:	2-8°C or -15 to -30°C (tolerates repeated freezing and thawing cycles)

Instructions for use

The Blocking Solution is ready-to use.

White-coloured fragments of casein can be found in the bottle after long-term storage at 4°C. They do not affect the quality and the blocking efficiency of *The Blocking Solution*. Please shake the buffer thoroughly before use.

After immobilization of capture antibody or target protein, *The Blocking Solution* is added without dilution to wells or membranes.

Incubate at room temperature for 1-4 hours or overnight (mostly one hour is quite enough). Duration of blocking can be minimized by shaking the plate at 600-900 rpm. Duration of blocking depends on characteristics of used surface and has to be tested individually.

After blocking, the surface has to be washed with *Washing Buffer* to make it useable for the next working steps.

After blocking, the coated antibodies or antigens can be stabilized by using *Liquid Plate Sealer*<sup>®</sup> (article number 160) for long-term storage at 4°C. After stabilization with *Liquid Plate Sealer*<sup>®</sup> a shelf life of several years can be observed with most antibodies.

Suitability of *The Blocking Solution* has to be tested in any assay, because of the different characteristics of assays and antibodies.

For further information please visit www.candor-bioscience.com.

Liquid Plate Sealer is a registered trade mark of CANDOR Bioscience GmbH.

